**PLANTOR**

**Kyiv-Kharkiv-Dnipro**

**2022**

**STOP THE BLEEDING - SAVE LIVES!**

**Hemostatic agent of local action with an antiseptic effect**

### Creation of "PLANTOR"

In connection with the spread of the war in Ukraine, an analysis of the data on the medical examination of the wounded and the causes of death in hostilities was carried out. One of the first is a massive blood loss. A group of Ukrainian scientists decided to create a hemostatic of local action (hemostop), which can be left to stay in the wound for a long time if evacuation of the victim is difficult or impossible. After careful analysis, the basis was particular substances of a wide range of actions are taken. The development process led to the creation of a unique tool, the main feature of which is the ability quickly stop bleeding of any type, including life-threatening bleeding, without causing further potential problems during surgery and in the rehabilitation period. The need for such means is very significant, especially in the conditions of hostilities, but the timely and effective stopping of bleeding saves life and is the guarantee of soldiers' return to the ranks. Different requirements for such means can prevent infections complication due to injury (antiseptic effect) and ease of use, including and in the order of self-help.

These are the main reason why the choice among hemostatics will be "PLANTOR".


"PLANTOR" was created and mass-produced by a team of scientists, technologists, and organizers under the general leadership of the director of the Official Representation in Ukraine of the International Nobel Information Center Dr. Sokolovsky S.I. based on the scientific development of the famous scientist Professor Shtrygol S.Yu. within the walls of the National Pharmaceutical University (Kharkiv) with the participation of the Academy of Technological Sciences of Ukraine (Kyiv) and LLC "Sokolovsky-Fortuna" (Dnipro).

### The principle of action of "PLANTOR"

Hemostatic "PLANTOR" is a tool of local action with an antimicrobial (antiseptic) effect, containing polysaccharides of vegetable origin, calcium salt, and miramistin. "PLANTOR" is safe even if left in the wound long time.

"PLANTOR" has a special sorption effect - the ability to fluid absorption, therefore, when in contact with blood absorbs plasma. Due to this, a blood clot is formed. At the same time, "PLANTOR" remains in the wound, does not combine with the blood, and does not enter the blood vessels. Moreover, its main component is natural polysaccharide carrageenan, which was selected as the most promising compared to other polysaccharides according to the results of laboratory studies. This means that "PLANTOR" exerts a powerful hemostatic effect. It does not affect the natural healing of the wound, does not slow down regeneration action, and does not change the chemical properties of tissues and the wound site. In the course of studies of "PLANTOR" and with its subsequent use, wounds were healed by primary tension, and normal scars were formed, which indicates the absence of a cytopathogenic effect.

When contacting those with blood "PLANTOR" particles begin to actively absorb liquid. In a matter of seconds, they swell and stick together into a single thick gelatinous mass that stops bleeding.

This clot sticks well to the edges of the wound, which literally covers the wound regardless of its type. The presence of an antiseptic component reliably prevents the development of wound infection. During further surgical treatment of the wound, the clot is easily and completely removed.

It is fundamentally important that "PLANTOR" does not interfere with the natural cascade mechanism of blood coagulation, and does not have a systemic effect - it only forms a clot from the blood with which it directly contacts the wound, due to the biophysical process of liquid absorption. "PLANTOR" does not contain procoagulant minerals or nanoparticles, such as smectite or kaolin, which remain in the human body until they are removed. Thus, it does not belong to means with a chemical or biochemical mechanism of action.

### Advantages of "PLANTOR":

"PLANTOR" has advantages that allow the manufacturer to position products as unique, and able to have an effect even in extreme conditions.

1. Able to stop bleeding for conditions states of hypothermia - a decrease in the internal temperature of the human body to a temperature lower than 35 C, when a large loss of blood quickly leads the victim to a state in which normal blood coagulation slows down and, in turn, it becomes difficult to stop bleeding.

"PLANTOR" is able to stop bleeding in hypothermic conditions because it acts independently of natural blood clotting factors.

2. Able to stop blood if it contains antiaggregants and anticoagulants (aspirin, warfarin, heparin, etc.) - substances that slow blood clotting. Despite their different principles of action, all of these drugs make it much more difficult to stop bleeding. "PLANTOR" is effective regardless of what medications the patient was taking (aspirin, heparin, and others).

3. It is a bioinert agent and does not generate heat. Even if "PLANTOR" particles remain in the wound, they do not harm the body, and do not cause inflammation. Unlike other hemostatics, when using "PLANTOR" the temperature in the wound increases by less than one degree Celsius, which makes it impossible to form a thermal or chemical burn.

4. "PLANTOR" in a complex with substances of synthetic origin (calcium gluconate) is capable of increasing the rate of onset of effective hemostasis by up to 10% in comparison with known means on the market of Ukraine.

5. "PLANTOR" is not toxic, does not cause allergic reactions compared to Celox, does not cause burns.

6. "PLANTOR" is cheaper than known hemostatic agents.

### Method of action and method of application of "PLANTOR"

"PLANTOR" is very easy to use: you just need to tear the bag, pour the contents on the wound and tightly bandage it. In the case of massive, severe bleeding (for example, arterial bleeding), it is necessary to combine the use of "PLANTOR" with mechanical methods of hemostasis (tourniquet, hemostatic tourniquet, etc.). If necessary, in case of heavy bleeding, it is recommended to use more than one dose.

Due to compactness, lightness and fast-acting use of "PLANTOR" is used in any emergency, critical situation without sterile conditions, as it provides an antiseptic effect.

### Purpose of "PLANTOR":

Stopping external bleeding of varying intensity:

- Everyone bleeding wounds.
- Lacerations, sores, and cuts.

### Features:

- Manufacturer: LLC "LekoPro" to order using the technology of LLC "Sokolovsky-Fortuna" with the scientific and organizational support of the Academy of Technological Sciences of Ukraine, the official representative office in Ukraine of the International Information Nobel Center, Dnipro Medical Institute of Traditional and Non-Traditional Medicine.
- Amount of drug in 1 sachet: 20 g.
- Significantly reduces blood loss.
- Valid even if anticoagulants are present in the blood, which prevent folding.
- Non-toxic, does not cause allergic reactions, or burns.
- Forms a blood clot (clot) in 30 seconds.
- Stops arterial bleeding.
- Made in Ukraine.
- Shelf life: 36 months from the date of manufacture.

After finishing after the expiration date, the powder can be used for up to 5 years, if the tightness of the package is not broken.

Ukraine is in dire need of hemostatic agents. This means "PLANTOR" saves lives in field conditions when the process of evacuation to hospital. Because these agents stop bleeding quickly in difficult conditions, they are extremely important. Therefore, its presence in every first-aid kit is necessary.

The hemostatic agent "PLANTOR" was created on volunteer principles and will be distributed free of charge Ukrainian military brigades until their full support.

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